

DATASHEET

The Deluxe Infant CRiSis Manikin features ECG monitoring, defibrillation, umbilical catheterization, training of i.v. and i.o. techniques as well as realistic CPR and airway management training.

The manikin is in size and details anatomically correct and ideal for practicing the pediatric advanced life support. Anatomical landmarks teach proper infant handling. Head, neck and jaw are fully articulated, allowing the jaw thrust maneuver. In addition, a manual brachial pulse is palpable.

Airway management

Train advanced airway management with the realistic anatomy of the mouth, tongue, oral pharynx, larynx, epiglottis, vocal cords, trachea and esophagus. The manikin has the capabilities of oral and digital intubation and allows the training of suction techniques.

Intravenous access on arm & leg

The injectable arms and legs have thin skin and the small lumen tubing simulates the size of an infant's vein. A pressurized system enables a realistic flashback.

Accessible are the cephalic and basilic veins on the arm, the great saphenous veins on the leg and the dorsal venous arches on the hand and foot.

Intraosseous infusion

With palpable landmarks including the patella, tibia and tibial tuberosity, and the pressurized system a realistic intraosseous infusion with aspiration of fluid can be trained. The skins and bones are replaceable.

ECG/Umbilical Cannulation Skin

This chest skin allows 4-lead ECG monitoring and the practice of umbilical cannulation, cutting and clamping. Fluids can be aspirated through the umbilicus.

Interactive ECG Simulator

With the Interactive ECG Simulator 17 pediatric heart rhythms can be recreated, as well as simulated pacing or cardioversion. It can be used with manual, semi-automatic or automatic defibrillators.

Deluxe Infant CRiSis™ Manikin



KEY BENEFITS

- Comprehensive PALS trainer
- Training of CPR and airway management
- 4-lead ECG monitoring
- With Interactive ECG Simulator: ECG simulation, defibrillation, pacing, cardioversion
- I.V. access, arm and leg
- I.O. access, leg
- Oral and digital intubation
- Suctioning capabilities
- Fully articulated head, neck and jaw
- Umbilical catheterization

Complete Infant CRiSis™ Manikin

The Complete Infant CRiSis Manikin includes the same features as the Deluxe Manikin, but is purchased without the Interactive ECG Simulator

ORDERING SPECIFICATIONS

Item no.	
LF03718U	Deluxe Infant CRiSis™ Manikin Fullbody Manikin with Interactive ECG Simulator, ECG/ Umbilical Cannulation Skin, IV bags, clamps, synthetic blood, winged infusion sets, syringes, 22-gauge needles, disposable umbilical cords and clamps, lubricant, Intraosseous bone kit with needle, lubricating jelly, towels, manual, carry case
LF03709U	Complete Infant CRiSis™ Manikin Fullbody Manikin with ECG/Umbilical Cannulation Skin, IV bags, clamps, synthetic blood, winged infusion sets, syringes, 22-gauge needles, disposable umbilical cords and clamps, lubricant, Intraosseous bone kit with needle, lubricating jelly, towels, manual, carry case
LF03670U	Interactive ECG Simulator



■ Ambu A/S
Baltorpbakken 13
2750 Ballerup
Denmark
T +45 72 25 20 00
F +45 72 25 20 50
ambu.com